

Creating Custom Prosthetic Arms Using Fusion 360's Solid and Mesh Modeling

MFG501520

Michael Peirone
Chief Executive Officer | @vichandproject

Kim Arklie
Mechanical Systems Designer | @vichandproject

# **About the Speakers**

- Michael Peirone
  - Chief Executive Officer, Victoria Hand Project (VHP)
  - B.Eng. Biomedical Engineering, University of Victoria, 2016
- Kim Arklie
  - Mechanical Systems Designer, Victoria Hand Project (VHP)
  - B.Eng. Mechanical Engineering, University of Victoria, 2021



# **Class Summary**

- The class will cover the following topics:
  - Victoria Hand Project's international work in the field of prosthetics
  - Issues with creating upper-arm (i.e above elbow) prosthetic devices
  - Combining solid bodies with 3D scan data to produce upper-arm prosthetic devices
  - Patient case study and workflow demonstration
  - How we implement our product workflows in partner clinics around the world

### **Learning Objectives**

- 1. Create custom healthcare devices by combining **3D scanning** and **CAD tools**
- 2. Implement custom product design workflows using the Fusion 360 API
- 3. Combine 3D scans and solid bodies using both Solid and Mesh Modeling in Fusion 360
- 4. Integrate Fusion 360 in clinical settings to improve productivity



### **International Work**





### **International Work**








### The Victoria Hand












3D-scan the patient's limb shape using **Autodesk Recap Pro** and take anatomical measurements




Casted Limb Impression



Autodesk Recap Pro




Generate the custom limb socket using **Autodesk Fusion 360** and **Meshmixer** 




Autodesk Meshmixer



3D-print the **custom** limb socket



Ultimaker Cura




# Trans-Humeral Prosthetics

**Background and Original Workflow** 

# **Development of Trans-Humeral Prosthetics**

### **Trans-Humeral Prosthetics Background**

- Help more amputees globally
- More complex scan data
- Design Constraints:
  - Good cosmetic appearance
  - Low weight
  - Appropriate fit and range of mobility
  - Functional elbow
  - Adaptable workflow



# **Original Trans-Humeral Socket Workflow**

### **Trans-Humeral Prosthetics Background**

- Used Autodesk Fusion 360 and Autodesk Meshmixer
- Various pre-sized lower sockets made in Fusion 360
- Patient's scan data smoothed and offset in Meshmixer to produce the upper socket
- Imported the pre-made lower socket into Meshmixer and combined with the upper socket




Fusion 360 Pre-Sized Lower Sockets

# **Original Trans-Humeral Socket Workflow**

### **Trans-Humeral Prosthetics Background**

- Used Autodesk Fusion 360 and Autodesk Meshmixer
- Various pre-sized lower sockets made in Fusion 360
- Patient's scan data smoothed and offset in Meshmixer to produce the upper socket
- Imported the pre-made lower socket into Meshmixer and combined with the upper socket



AutoDesk MeshMixer Socket Workflow Video: https://www.youtube.com/watch?v=Fo0p7GjkJxU

### **Issues with the Original Trans-Humeral Socket Workflow**


### **Trans-Humeral Prosthetics Background**

- 1. Not a smooth transition between lower and upper portions of the socket
- 2. Cannot resize the socket easily due to the cross-platform workflow (i.e. export from Fusion 360, import into Meshmixer)
- 3. Too many user inputs and commands



Poor Quality Socket Transition






# Trans-Humeral Prosthetics

Optimization of Workflow

### **Workflow Overview**


- Utilizes the Mesh Workspace in Autodesk Fusion 360
- Pre-made lower socket was redesigned to be parametrically dimensioned, allowing for more customization
- Limb impression is imported into Fusion 360 as a mesh and kept as a mesh body to be edited



### **Optimization of Trans-Humeral Socket Workflow**

# 1. Create the Lower Socket


Size the lower socket by inputting the patient's recorded measurements into the solid body design parameters.



**Optimization of Trans-Humeral Socket Workflow** 

# 2. Import and Align the Limb Impression


Insert the limb impression into Fusion 360 as a mesh body. Upon insertion, align the limb impression to the lower socket solid body.



### **Optimization of Trans-Humeral Socket Workflow**

# 3. Create the Upper Socket


Duplicate, trim, repair, remesh, smooth, and offset the limb impression 5mm outwards to create the basis of the upper socket.



### **Optimization of Trans-Humeral Socket Workflow**

# 4. Tessellate the Lower Socket


Confirm the alignment of the limb impression, and lower and upper socket. Then tessellate the lower socket into a mesh body.



### **Optimization of Trans-Humeral Socket Workflow**

# 5. Boolean Union the Lower and Upper Socket


Perform a Boolean Union operation to combine the previously prepared lower socket and 5mm offset of the upper socket.



### **Optimization of Trans-Humeral Socket Workflow**

# 6. Clean Up Socket Transition


Combine the face groups on the main socket body. Erase and fill, remesh, and smooth the socket to achieve a seamless transition between the two socket components (i.e. lower and upper).



### **Optimization of Trans-Humeral Socket Workflow**

# 7. Boolean Difference the Limb Impression and Combined Socket

Perform a Boolean Difference operation to cut the limb impression from the combined lower and upper socket. The socket is now complete and ready to be prepared for 3D printing.




# **Advantages**

- Higher precision model
- Improved socket fit for better usability of the prosthetic device



# **Advantages**




- Uniform thickness of the upper socket
- Improved socket appearance
- Integration with the Fusion 360 API

# **Advantages**

- Semi-automate alignment process
- More easily modify and realign components throughout workflow
- Simplification of workflow







- 7.0 magnitude earthquake occurred in Haiti on January 12, 2010 [3]
- Hundreds of thousands were killed or injured in the disaster [3]




- John was injured in the earthquake when his house collapsed on top of him
- Never received a prosthetic arm
- Recently met with VHP in-country partner clinician in Haiti to be fit with a device






- In-country clinician completes end-toend workflow independently
- Plaster limb impression and patient measurements taken






- Limb impression is digitized using photogrammetry methods
- Optimized trans-humeral upper socket workflow is followed to create the upper socket
- Clinician applied their knowledge and expertise to create a proper fitting socket



#### **Haiti Trans-Humeral Patient**

Socket is exported from Fusion 360 and 3D printed on an Ultimaker 3D printer


 Victoria Hand and additional trans-humeral components are 3D printed by the local 3D Print Technician



- Trans-humeral system is assembled by the local 3D Print Technician
- Clinician fits patient with the upper socket and prosthetic system
- In-country partner team provides feedback to VHP team in Canada







## Workflow Implementation

**VHP Partner Clinics** 

## **Workflow Deployment Model**

**Workflow Implementation in VHP Partner Clinics** 

- Workflow provided to VHP's partners
- Integrated into VHP software
- Easier to produce custom devices




Victoria Hand Project. Egypt, 2019

# Combining Advanced Technologies with Traditional Workflows

#### **Workflow Implementation in VHP Partner Clinics**

- Combine CAD tools and additive manufacturing with traditional prosthetic fitting methodologies
- More efficiently create devices at a lower cost
- Allow clinicians to focus more time on patient care
- Produce better custom fit devices



Victoria Hand Project Software Demo

Video: https://www.youtube.com/watch?v=tC7kadLzsJU

# Combining Advanced Technologies with Traditional Workflows

#### **Workflow Implementation in VHP Partner Clinics**

- Combine CAD tools and additive manufacturing with traditional prosthetic fitting methodologies
- More efficiently create devices at a lower cost
- Allow clinicians to focus more time on patient care
- Produce better custom fit devices.



### **Future Work**

#### **Workflow Implementation in VHP Partner Clinics**

- Refining the development of the transhumeral socket workflow
- Automating the workflow further with VHP software and the Fusion 360 API to make the process of creating the socket easier for the clinician
- Realistic color information for increased accuracy when creating the socket



## **Future Work**

#### **Workflow Implementation in VHP Partner Clinics**


- Refining the development of the transhumeral socket workflow
- Automating the workflow further with VHP software and the Fusion 360 API to make the process of creating the socket easier for the clinician
- Realistic color information for increased accuracy when creating the socket



Socket Shape Sketched on Impression


## **To Summarize**

- Custom 3D printed prosthetic arms
- Previous trans-humeral socket workflows were too tedious and resulted in sockets with a poor appearance
- Use Fusion 360 Mesh Workspace to easily create custom prosthetic sockets
- Easier workflow resulting in a socket with a better cosmetic appearance



## **To Summarize**

- Custom 3D printed prosthetic arms
- Previous trans-humeral socket workflows were too tedious and resulted in sockets with a poor appearance
- Use Fusion 360 Mesh Workspace to easily create custom prosthetic sockets
- Easier workflow resulting in a socket with a better cosmetic appearance



## **Our Team**



Kelly Knights, B. Eng., EIT

**Biomedical Systems Designer** 



Jacqui Moreland, B. Eng.

Biomedical Software Developer



Dr. Nick Dechev, Ph.D., P. Eng.

**Chief Technical Officer** 





**Prosthetics Within Reach** 



Victoria Hand Project



@VicHandProject



@VicHandProject



Victoria Hand Project



Victoria Hand Project



www.victoriahandproject.com

### References

[1] B. Phillips, G. Zingalis, S. Ritter, and K. Mehta, "A review of current upper-limb prostheses for resource constrained settings," 2015 IEEE Global Humanitarian Technology Conference (GHTC), 2015.

[2] M. LeBlanc, "Give Hope - Give a Hand - The LN-4 Prosthetic Hand," November 2008. [Online]. Available: <a href="https://web.stanford.edu/class/engr110/2011/LeBlanc-03a.pdf">https://web.stanford.edu/class/engr110/2011/LeBlanc-03a.pdf</a>.

[3] "2010 Haiti earthquake," History.com, 18-Jul-2011. [Online]. Available: https://www.history.com/this-day-in-history/massive-earthquake-strikes-haiti. [Accessed: 19-Aug-2022].

